THERMAL TRENDS

6

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- L75 PT1600
- STA PT1000
- STA PT1600
- DSC PT1600
- TMA PT100
- Thermo KINETICS

Dilatometer L75 PT1600

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
Operation Modes

Standard	Included	Included	Included
Stress/Strain		Included	Optional
Creep		Included	Optional
Stress Relaxation		Included	Optional
Dynamic TMA		Included	Optional

Optional

Dta feature	Optional	Optional	Optional
-------------	----------	----------	----------

ADVANCED TOOLS FOR THERMAL AGEING AND SAFETY ANALYSIS

B. Roduit⁽¹⁾, Ch. Borgeat⁽¹⁾, B. Berger⁽²⁾, P. Folly⁽²⁾, J.N. Aebischer⁽³⁾, H. Fierz⁽⁴⁾ and F. Stoessel⁽⁴⁻⁵⁾

⁽¹⁾ Advanced Kinetics and Technology Solutions (AKTS AG) TECHNO-Pôle, CH-3960 Siders - <http://www.akts.com>

⁽²⁾ armasuisse, Science and Technology Centre, CH-3602 Thun - <http://www.armasuisse.ch>

⁽³⁾ University of Applied Sciences of Western Switzerland, CH-1705 Fribourg - <http://www.eif.ch>

⁽⁴⁾ Swiss Institute for the Promotion of Safety & Security (SWISSI), CH-4002 Basel - <http://www.swissi.ch>

⁽⁵⁾ Swiss Federal Institute of Technology (EPFL), Institute of Process Sciences, CH-1015 Lausanne - <http://www.epfl.ch>

INTRODUCTION

Generally, all energetic materials evolve heat during decomposition. Processing, design, quality control, and operational applications all require an understanding of thermal hazards and an ability to predict safety limits and decomposition process in extended temperature ranges [1,2].

Several methods have been presented for predictions of the reaction progress of exothermic reactions under heat accumulation conditions [3]. However, because decomposition reactions usually have a multi-step nature, the accurate determination of the kinetic characteristics strongly influences the ability to correctly describe the progress of the reaction [4,5]. The use of simplified and conservative kinetic models for the assessment of runaway reactions leads to economic drawbacks, since they result in unnecessary large safety margins.

Applying the results obtained by :

- DTA (Differential Thermal Analysis)
- DSC (Differential Scanning Calorimetry)
- TG (Thermogravimetry)
- EGA (Evolved Gas Analysis MS or FTIR)

advanced numerical techniques such as AKTS-Thermokinetics and AKTS-Thermal Safety Software enable prediction of the reaction progress of materials in broad temperature range. In fact, numerical simulations are used to replace experiments in situations, which are not directly accessible to the experiment for timing or safety reasons.

ANALYSIS PROCESS

As a general rule, solid state reactions demonstrate profound multi-step characteristics as presented in figure 1. The assumption that the decomposition of an energetic material will obey a simple rate law is not often true. Using for example DSC data, the analysis process requires determination of the kinetic characteristics of the reaction.


Figure 1: Analysis process. Normalized DSC-signals of a nitrocellulose based compound as a function of the temperature for a reaction involving at least two exothermic events. Experimental data are represented as symbols, solid lines represent the calculated signals. The values of the heating rate in °C/min are marked on the curves.


Figure 2: Friedman analysis of the examined material.

Friedman analysis, based on the Arrhenius equation, applies the logarithm of the conversion rate dx/dt as a function of the reciprocal temperature at different degrees of the conversion.

$$\ln \frac{dx}{dt} = \ln(A) - \frac{E}{RT_{i,j}} + \ln(f(x_i))$$

with i : index of conversion, j : index of heating rate, $f(x)$ the function dependent on the decomposition mechanism. As $f(x)$ is constant at each conversion degree x_i , the method is so-called 'isoconversional'.


Figure 3: Activation energy as a function of the reaction progress for decomposition of the high energetic material (DSC closed crucibles).

The accurate determination of the kinetic parameters and optimization of the baseline which enable the correct fit of the experimental data is a prerequisite for prediction of the reaction progress under any new temperature profile (see figure 4).


Figure 4: Prediction of the reaction extent (DSC, normalized signals) and confidence interval of a nitrocellulose based substance as a function of time under isothermal conditions ($T = 120^\circ\text{C}$). These values indicate that there is a 95% probability that the reaction progress after 10 days exposure at 120°C is greater than 77 and lower than 93%. These values are in good agreement with a subsequent measurement under isothermal conditions (82% decomposition at 120°C after 10 days).

TMR_{ad} SAFETY MARGINS & HEAT ACCUMULATION CONDITIONS

The calculated kinetic parameters can be subsequently employed to predict the reaction progress of the investigated samples under any given temperature mode. For chemical process safety adiabatic conditions are used for the prediction of the Time to Maximum Rate under adiabatic conditions (TMR_{ad}).


Figure 5: Adiabatic runaway curves showing the confidence interval for the prediction ($T_{\text{beg}}=100^\circ\text{C}$ and $\Delta T_{\text{ad}}=\Delta H/c_p=2578\pm 173^\circ\text{C}$).

The second field of application for numerical simulation techniques in process safety is the solution of partial differential equations as they are encountered in heat conduction problems. These problems arise when heat accumulation situations are to be analyzed [6].

$$\frac{dT}{dt} = \frac{\lambda}{c_p \rho} \left(\frac{\partial^2 T}{\partial x^2} + \frac{\partial^2 T}{\partial y^2} \right) \quad (\text{cartesian}) \quad \frac{dT}{dt} = \frac{\lambda}{c_p \rho} \left(\frac{\partial^2 T}{\partial r^2} + \frac{1}{r} \frac{\partial T}{\partial r} \right) \quad (\text{cylindrical})$$

Applications of Finite Element Methods (FEM) and accurate kinetic description enable the determination of the effect of scale, geometry, heat transfer (isolation), thermal conductivity and ambient temperature on the heat accumulation conditions. In fact, the assumption that it is safe to handle an energetic material at any temperature below the first appearance of an exothermic signal on the DSC curve can be often false. The highest safe temperature for handling any energetic material depends on several factors such as its size, shape, and previous thermal history. Due to insufficient thermal convection and limited thermal conductivity, a progressive temperature increase in the sample can easily take place, resulting in a thermal explosion. Safe operating conditions with tailored safety margins can be defined using numerical simulation. Presented examples illustrate a slow cook-off experiment and the calculations performed on pressed PBX.


Figure 6: Slow cook-off experiments of pressed PBX based on 94% RDX (A) and simulation (B). As presented in the inset, the predicted temperature of explosion was 181°C . It is in good agreement with the slow cook-off experiments (179°C).


Figure 7: 3D calculations representing the different stages of the slow cook-off performed on pressed PBX.

CONCLUSIONS

Employing effective and efficient mathematical modeling, advanced kinetic analysis enables the calculations of the progress of decomposition reactions under temperature conditions different from those at which the original examinations were carried out. Applications of Finite Element Methods (FEM) and accurate kinetic description enable determination of the effect of scale, geometry, heat transfer, thermal conductivity and ambient temperature on the heat accumulation conditions.

REFERENCES

- [1] Stoessel F, Steinbach J, Eberz A: Plant and process safety, exothermic and pressure inducing chemical reactions. In: Ullmann's encyclopedia of industrial chemistry. Weise E (Eds), VCH, Weinheim (1995):343-354.
- [2] Keller A, Stark D, Fierz H, Heinze E, Hungerbuehler K: Estimation of the TMR using dynamic DSC experiments. Journal of Loss Prevention in the Process Industries (1997) 10(1):31-41.
- [3] D.A. Frank-Kamenetskii, Diffusion and Heat Transfer in Chemical Kinetics, Plenum Press, New York, London, 1969.
- [4] B. Roduit, Prediction of the progress of solid state reaction under different temperature modes, Thermochim. Acta, 388 (2002) 377.
- [5] B. Roduit, Computational aspects of kinetic analysis. Part E: The ICTAC Kinetics Project - Numerical techniques and kinetics of solid state processes, Thermochim. Acta 355, 171-180, 2000.
- [6] Gygax, R., Thermal Process Safety, Data Assessment, criteria, measures, ed. ESCIS. Vol. 8, 1993, Lucerne: ESCIS.
- [7] <http://www.akts.com> (AKTS-Thermokinetics and AKTS-Thermal Safety Software)

LINSEIS

LINSEIS GmbH	LINSEIS Inc.
Vielitzerstr. 43	20 Washington Road
95100 Selb	P.O.Box 666
Germany	Princeton-Jct. NJ 08550
Tel.: (+49) 9287-880 - 0	Tel.: (609) 799-6282
Fax: (+49) 9287-70488	Fax: (609) 799-7739
E-mail: info@linseis.de	E-mail: info@linseis.com

Products: DIL, TG, STA, DSC, HDSC, DTA, TMA, MS/FTIR, Laser Flash
Services: Service Lab, Calibration Service

www.linseis.com