


DIL Q PLATINUM



thermal analysis
with **out** limits



LINSEIS

TA



This dilatometer is a very unique instrument, which was developed for customers in the ceramics automotive industry.

These customers have to make a lot of quality control measurements in order to ensure ISO9000 standards. Especially companies manufacturing catalytic converter base ceramics for car exhaust systems are very interested.

The Quattro Dilatometer is built up with four separate dilatometer measuring sensors, which can measure simultaneously either four separate samples at one time, or three separate samples against a NIST reference.

That means, the productivity of the Quattro Dilatometer is three times as high, if compared to the normally used dual push rod dilatometers. With a dual push rod dilatometer only one sample can be measured against the standard at each time.

In order to further increase the productivity, this dilatometer is built with two separate furnaces. Whenever the first run of four samples is ready, a new, cool furnace is ready to immediately start the next run.

There is another feature included in this Quattro Dilatometer, which is called an automatic furnace lift mechanism. This feature automatically lifts the furnace at the end of each measurement without any operator interaction.

This way the measuring system is already cool when the operator comes to change the samples. This again increases the productivity of the system.

The 32 bit software was further developed to be able to program the parameters for four samples on one screen. For each sample different sample length, file name etc can be stored, and are available for later evaluation.

The temperature programming is done through an extensive program part, for several stages, dwell times, heat up speeds etc. A software macro is supplied with automatic evaluation, to ensure that ready measuring data are available without time delay.

We are proud to say that as far as we know Linseis is the only supplier worldwide, that can offer this highly specialized dilatometer. This is another example of our powerful line of thermal analysis instruments.

www.linseis.com

International

LINSEIS Messgeräte GmbH
Vielitzerstrasse 43
95100 Selb / Germany
Tel.: 09287 / 880 - 0
Fax.: 09287 / 70488
E-mail: linseis@t-online.de



USA

Linseis Inc.
20 Washington Road
P.O.Box 666
Princeton-Jct.NJ 08550
Tel.: (609) 799 - 6282
Fax.: (609) 799 - 7739
E-mail: info@linseis.com



L I N S E I S